





Marking duplicates

Removing non-independent observations

Data Pre-processing for Variant Discovery




Mark duplicates to mitigate duplication artifacts



Duplicates = **non-independent measurements** of a sequence fragment


-> Must be removed to assess support for alleles correctly



✗ = sequencing error propagated in duplicates

Why mark duplicates?



- **Non-independent measurements** of sequence
 - Sampled from single template of DNA
 - Violates independence assumptions made in variant calling
- Errors in sample/library prep are propagated to *all* the duplicates
- “Best” copy – mitigates the effects of errors



✗ = library prep error propagated in duplicates

Where does the duplication come from?


- **LIBRARY DUPLICATES**
 - Increases with PCR cycles
- **OPTICAL DUPLICATES**
 - Are nearby clusters on a flow cell lane




Optical and PCR duplication events arise at different rates as a sequencing experiment proceeds




PCR duplicates



Optical duplicates



Duplicates are flagged the same but can be tagged differently (DT)



How do we identify duplicate reads?



- Dupes might come from the same input DNA template, so we will assume that reads will have same start position on reference
 - “Where was the first base that was sequenced?”
 - For paired-end (PE) reads, same start for both ends
- Identify duplicate sets, then choose representative read based on base quality scores and other criteria

But there's a catch (or two)...



- BWA sometimes “clips” bases from the ends of the alignment (when the alignment there is poor)
 - Need to use SAM flags + CIGAR string to determine the unclipped 5' end
- Fragments mapped to the reverse strand are specified by their 3' position, instead of 5'

Identify duplicates using orientation + “unclipped” 5’ position

Pos	1	2	3	4	5	6	7	8	9
Ref	T	A	G	C	C	G	A	T	C
r1	<u>T</u>	A	G	C	C	G	A		
r2	T	A	G	C	C	G	<u>A</u>		
r3	<u>T</u>	A	—	C	CAG	A			
r4	T	A	G	C	C	H	<u>H</u>		
r5	<u>T</u>	A	G	C	C	G	A	T	C
r6	<u>S</u>	S	G	C	C	G	A		
r7			<u>G</u>	C	C	G	A		

Blue maps to forward strand

Red maps to reverse strand

Grey bases are clipped

Underlined is the expected 5’ start of the read, given the mapping

What are the duplicate sets?

Identify duplicates using orientation + “unclipped” 5’ position

Pos	1	2	3	4	5	6	7	8	9
Ref	T	A	G	C	C	G	A	T	C
r1	<u>T</u>	A	G	C	C	G	A		
r2	T	A	G	C	C	G	<u>A</u>		
r3	<u>T</u>	A	—	C	CAG	A			
r4	T	A	G	C	C	H	<u>H</u>		
r5	<u>T</u>	A	G	C	C	G	A	T	C
r6	<u>S</u>	S	G	C	C	G	A		
r7			<u>G</u>	C	C	G	A		

Blue maps to forward strand

Red maps to reverse strand

Grey bases are clipped

Underlined is the expected 5’ start of the read, given the mapping

So...what are the duplicate sets?

☞ r1, r3, r5, r6 (start at position 1)

Identify duplicates using orientation + “unclipped” 5’ position

Pos	1	2	3	4	5	6	7	8	9
Ref	T	A	G	C	C	G	A	T	C
r1	<u>T</u>	A	G	C	C	G	A		
r2	T	A	G	C	C	G	<u>A</u>		
r3	<u>T</u>	A	—	C	C	A	G	A	
r4	T	A	G	C	C	H	<u>H</u>		
r5	<u>T</u>	A	G	C	C	G	A	T	C
r6	<u>S</u>	S	G	C	C	G	A		
r7			<u>G</u>	C	C	G	A		

Blue maps to forward strand

Orange maps to reverse strand

Grey bases are clipped

Underlined is the expected 5’ start of the read, given the mapping

So...what are the duplicate sets?

☞ r1, r3, r5, r6 (start at position 1)

☞ r2, r4 (start at position 7)

Identify duplicates using orientation + “unclipped” 5’ position

Pos	1	2	3	4	5	6	7	8	9
Ref	T	A	G	C	C	G	A	T	C
r1	<u>T</u>	A	G	C	C	G	A		
r2	T	A	G	C	C	G	<u>A</u>		
r3	<u>T</u>	A	—	C	CAG	A			
r4	T	A	G	C	C	H	<u>H</u>		
r5	<u>T</u>	A	G	C	C	G	A	T	C
r6	<u>S</u>	S	G	C	C	G	A		
r7			<u>G</u>	C	C	G	A		

Blue maps to forward strand

Orange maps to reverse strand

Grey bases are clipped

Underlined is the expected 5’ start of the read, given the mapping

So...what are the duplicate sets?


☞ r1, r3, r5, r6 (start at position 1)

☞ r2, r4 (start at position 7)


☞ r7 (starts at position 3)

So now we have mapped, sorted, and *deduped* reads

Showing duplicate reads



Hiding duplicate reads



- Duplicate status is indicated in SAM flag
- Duplicates are **not removed**, just tagged (unless you request removal)
- Downstream tools can read the tag and choose to ignore those reads
- Most GATK tools ignore duplicates by default

Use cases where you may *NOT* want to mark duplicates




- Amplicon sequencing (or any other LC with non-random starts)
 - All reads start at same position by design
 - If have UMI can use `UMIAwareMarkDuplicatesWithMateCigar`
- RNAseq (allele-specific) expression analysis
(alternatively, `ASEReadCounter` can disable `DuplicateFilter`)

Add-on: Predicting the complexity of a sequencing experiment

Complexity analysis depends on:

- Estimated library size
- Return on Investment (ROI) calculations



Estimation of library size and duplication in Picard

Mathematical Notes on SAMtools Algorithms

Heng Li

October 12, 2010

Duplicate Rate

1.1 Amplicon duplicates

Let N be the number of distinct segments (or seeds) before the amplification and M be the total number of amplicons in the library. For seed i ($i = 1, \dots, N$), let k_i be the number of amplicons in the library and k_i is drawn from Poisson distribution $\text{Po}(\lambda)$. When N is sufficiently large, we have:

$$M = \sum_{i=1}^N k_i = N \sum_{k=0}^{\infty} k p_k = N \lambda$$

where $p_k = e^{-\lambda} \lambda^k / k!$.



Estimated fraction of
duplicates

$$d \simeq 1 - \frac{N}{m} \left(1 - e^{-m/N} \right)$$

Assumptions

- all reads are drawn from the same Poisson distribution $\text{Po}(\lambda)$
- the occurrence of duplication events depends on underlying concentration of inserts in the library

Active research to improve library size estimation



Coverage Bias and Sensitivity of Variant Calling for Four Whole-genome Sequencing Technologies

Nora Rieber^{1*}, Marc Zapatka^{2*}, Bärbel Lasitschka³, David Jones⁴, Paul Northcott⁵, Barbara Hutter¹, Natalie Jäger¹, Marcel Kool⁶, Michael Taylor^{5,6}, Peter Lichter², Stefan Pfister^{4,7}, Stephan Wolf³, Benedikt Brors¹, Roland Eils^{1,8*}

- Rate of duplication varies with insert size length
- Duplications rates also likely vary with GC content

Data Pre-processing for Variant Discovery

