




Somatic SNV and Indel Discovery


Basic operation and algorithm of
Mutect2 and related tools

Somatic SNV & Indel discovery workflow



Logic of the Tumor-Normal workflow

Comparison to matched normal -> subtraction of germline background



Tumor-only analysis



- It is possible to run the workflow without a matched normal in “tumor-only mode” (normally used for PON creation)
- MUST have a good PON to eliminate common germline variation
- Will still require extra filtering (not described here)





Panel of Normals for SNVs & Indels




- VCF of calls made from a set of unrelated “normal” samples
- Main purpose:
Eliminate common/recurring technical artifacts
-> should use normals made using the same data generation techniques
(eg same capture kit for exomes, same sequencing platform etc)
- Secondary purpose: also eliminates germline variants not called in the matched normal (or approximates the normal if none is available)

False positives from artifacts and germline variation

Somatic point mutations occur ~ 1 / Mbp

	ARTIFACT	GERMLINE EVENT
TUMOR		
NORMAL		
At risk	Every base	~ 1667 germline variants / Mbp
Source	<ul style="list-style-type: none"> Misread bases Sequencing artifacts Misaligned reads 	<ul style="list-style-type: none"> Low coverage in NORMAL
Solutions	<i>filters, Panel of Normals (PoN)</i>	<i>dbSNP, ExAC, COSMIC, PoN</i> gnomAD



Step 1

CALL VARIANTS WITH MUTECT2


Mutect2 is based on HaplotypeCaller

Skip :

- Sites in *PoN*
- Sites with high fraction alt alleles in normal

Allele-specific calling:

- Distinguishes alleles in the *germline population frequency* resource and uses AF in calculating probability variant exists in normal *and* tumor



Somatic genotypes inferred from PairHMM likelihoods

$$\mathcal{L}(G_{\text{ref}}|R) = \prod_j \mathcal{L}(G_{\text{ref}}|R_j) \quad \mathcal{L}(G_i|R) = \prod_j [\mathcal{L}(G_i|R_j)f_{\text{alt}} + \mathcal{L}(G_{\text{ref}}|R_j)(1 - f_{\text{alt}})]$$

Likelihood of reference
genotype given all reads

Likelihood of variant
genotype i given all reads

Likelihoods of variant/reference
alleles given read j

$$\text{LR}_i = \log \mathcal{L}(G_i|R) - \log \mathcal{L}(G_{\text{ref}}|R)$$

Log-likelihood ratio for
genotype i


- No explicit ploidy assumptions (unlike HaplotypeCaller)
 - somatic genotype likelihoods weighted by variant allele fraction
- Statistical threshold for somatic call uses log-likelihood ratios
 - ≥ 5.3 in favor of the variant somatic genotype
 - Also filter based on the likelihood of the allele in the Normal

Converting likelihoods to probabilities uses AF (f)

- If variant is in gnomAD:
Use the allele frequency f in gnomAD
- If variant is *not* in gnomAD:
*Set $f = 1/(2 * \text{\#samples} + 2)$, which is a reasonable guess for the allele frequency of the variant*
- If we don't have gnomAD:
Default to $f = 0.001$

Genotyping in GATK4 Mutect2

- GATK4 Mutect2 models the allele fractions and allele assignment to each read as latent variables f and z
- Choose the allele set A that maximizes model evidence
- If log odds > 3.0 (by default) then emit variant
- At low coverage sites, the Bayesian approach in GATK4 performs better than the prior frequentist approach in GATK3



$$f \sim \text{Dirichlet}(\alpha)$$

$$z|f \sim \text{Categorical}(f)$$

$$p(r|z_{ra}) = l_{ra}$$

$$l_{ra} \equiv P(\text{read } r | \text{allele } a)$$

from PairHMM


$$\log \frac{p(\mathbb{R} | A_{alt})}{p(\mathbb{R} | A_{ref})} > \delta = 3.0$$

then emit variant

Case Study: 120 base deletion




Tumor: BWA alignment



Clear evidence of some sort of event is present, but impossible for a traditional pileup-based caller to recover

Normal: BWA alignment




Event would also not be caught with discordant read pair caller, since insert sizes of supporting reads are within normal range

Case Study: 120 base deletion



Tumor: M2 realignment




MuTect2
reassembly
recovers the
120 base
deletion
haplotype

Normal: M2 realignment




It also discerns
reads that are
unambiguously
phased into the
WT haplotype,
and a haplotype
with insufficient
likelihood.



Step 2

FILTER RAW VARIANT CALLS

Filtering is based on annotations + contamination estimate



ANNOTATION	INFO field annotations
Coverage	DP
DepthPerAlleleBySample	AD
TandemRepeat	STR
OxoGReadCounts	F1R2, F2R1
ReadPosition	MPOS
BaseQuality	MBQ
MappingQuality	MMQ
FragmentLength	MFRL
StrandArtifact	SA_POST_PROB, SA_MAP_AF

** Not a comprehensive list*

FilterMutectCalls filters for multiple criteria

FILTER	Description
artifact_in_normal	artifact_in_normal
base_quality	alt median base quality
clustered_events	Clustered events observed in the tumor
contamination	contamination
duplicate_evidence	evidence for alt allele is overrepresented by apparent duplicates
fragment_length	abs(ref - alt) median fragment length
germline_risk	Evidence indicates this site is germline, not somatic
mapping_quality	ref - alt median mapping quality
multiallelic	Site filtered because too many alt alleles pass tumor LOD
orientation_bias	Orientation bias (in one of the specified artifact mode(s) or complement) seen in one or more samples.
panel_of_normals	Blacklisted site in panel of normals
read_position	median distance of alt variants from end of reads
str_contraction	Site filtered due to contraction of short tandem repeat region
strand_artifact	Evidence for alt allele comes from one read direction only
t_lod	Tumor does not meet likelihood threshold

Additional filters for sequence context artifacts

FILTER	Description
artifact_in_normal	artifact_in_normal
base_quality	alt median base quality
clustered_events	Clustered events observed in the tumor
contamination	contamination
duplicate_evidence	evidence for alt allele is overrepresented by apparent duplicates
fragment_length	abs(ref - alt) median fragment length
germline_risk	Evidence indicates this site is germline, not somatic
mapping_quality	ref - alt median mapping quality
multiallelic	Site filtered because too many alt alleles pass tumor LOD
orientation_bias	Orientation bias (in one of the specified artifact mode(s) or complement) seen in one or more samples.
panel_of_normals	Blacklisted site in panel of normals
read_position	median distance of alt variants from end of reads
str_contraction	Site filtered due to contraction of short tandem repeat region
strand_artifact	Evidence for alt allele comes from one read direction only
t_lod	Tumor does not meet likelihood threshold

→ FilterByOrientationBias

E.g. likely OxoG G→T transversions

Mutect2 command and main options

Base command for PoN creation and tumor-only analysis:

```
gatk Mutect2 \  
  -R ref_fasta.fa \  
  -I sample.bam \  
  -tumor sample_name \  
  -L intervals.list \  
  -O sample.vcf.gz
```

For matched-normal tumor calling add:

```
-I normal.bam \  
-normal normal_sample_name \  
-bamout bamout.bam \
```

Reassembled BAM now recommended

To specify a germline AF resource:

```
--germline_resource af-only-gnomad.vcf.gz \  
--af_of_alleles_not_in_resource 0.0000025 \
```

Germline resource must have allele-specific frequencies; af for not in gnomAD exomes

To specify a PoN:

```
--normal_panel pon.vcf.gz \
```

Filtering commands and main options

Filter M2 calls for multiple contexts:

```
gatk FilterMutectCalls \  
  -V tumor_matched_m2_snvs_indels.vcf.gz \  
  -contaminationTable contamination.table \  
  -O tumor_matched_m2_oncefiltered.vcf.gz
```

*Output of CalculateContamination;
FilterMutectCalls uses the first row listing BAM
file-level contamination*

Afterwards, optionally filter by orientation bias:

```
gatk FilterByOrientationBias \  
  -V tumor_matched_m2_oncefiltered.vcf.gz \  
  --artifactModes 'G/T' \  
  -P tumor.preadapter_detail_metrics \  
  -O tumor_oxog_twicefiltered.vcf.gz
```

*Requires pre-adapter detailed metrics calculated
by Picard CollectSequencingArtifactMetrics.*

Somatic SNV & Indel discovery workflow

